

# CSSD Checklist for NABH Compliance

## ● General Infrastructure

- Separate entry and exit for sterile and unsterile items
- Ventilation system to prevent cross-contamination
- Sufficient lighting and space for workflow
- Clearly defined zones: Decontamination, Packing, Sterilization, Storage

## ● Sterilization Equipment

- Autoclave with automatic controls and temperature monitoring
- Biological and chemical indicators used for each sterilization cycle
- Regular maintenance and calibration records of equipment
- Logbooks for each cycle and load details maintained

## ● Instrument Washing & Preparation

- Designated instrument washing sink with Sliding Instrument Trays
- Use of enzymatic solutions and soft brushes
- Proper drying, inspection and preparation of trolleys

## ● Trolleys and Transportation

- Dedicated CSSD trolleys for sterile and unsterile material
- Clearly labeled and easy to clean surfaces
- Segregation of clean and contaminated items during transport

## ● Storage and Organization

- Sterile storage area with temperature and humidity control
- Sterile items stored above the floor level and away from walls
- Use of Lapro Instrument Sterilization Box cum Organizer
- FIFO (First In, First Out) method followed for stock rotation

\*\*The above checklist is for reference purposes only.